CELLULOSE LACTOSE

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| Sieve data – co-processed lactose | | |
|-----------------------------------|--------------|-------------------------------------|
| | Lactose type | Cellactose® 80 specified/typical |
| Particle size distribution | < 32 µm | NMT 20%/7% |
| Method: | < 160 µm | 35 – 65%/54% |
| Air jet sieving | < 250 µm | NLT 80%/93% |

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Functional related characteristics

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
Compactability

High-dose Vitamin C formulation with Cellactose® 80

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

Adherence capacity

Content uniformity – Cellactose® 80 vs. various excipients

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Figure 14: MEGGLEのCellactose® 80の包装、品質保証期間。

Literature

- [1] Meeus, L. (2011). Direct Compression versus Granulation. *Pharmaceutical Technology*, 23(3).
- [2] Kristensen, H. G., & Schaefer, T. (1987). Granulation: A Review on Pharmaceutical Wet-Granulation. *Drug Development and Industrial Pharmacy*, 13(4–5), 803–872.
- [3] Mîinea, L. A., Mehta, R., Kallam, M., Farina, J. A., & Deorkar, N. (2011). Evaluation and Characteristics of a New Direct Compression Performance Excipient, 35(3).
- [4] Schmidt and Rubensdörfer (1994). Evaluation of Ludipress as a “Multipurpose Excipient” for DC Part I: Powder Characteristics and Tableting Properties, *Drug dev. ind. Pharm.* 20(18), 2899–2925

MEGGLE App:

MEGGLE Consultant

MEGGLE Group Wasserburg
BG Excipients & Technology
Megglestrasse 6–12
83512 Wasserburg
Germany

Phone +49 8071 73 476
Fax +49 8071 73 320
service.pharma@meggle.de
www.meggle-pharma.com

MEGGLE warrants that its products conform to MEGGLE's written specification and makes no other expressed or implied warranties or representations. For any specific usage, the determination of suitability of use or application of MEGGLE products is the sole responsibility of the user. The determination of the use, application, and compliance of this product with regard to any national, regional, or local laws and/or regulations is the sole responsibility of the user, and MEGGLE makes no representation with regards to same. Nothing herein shall be construed as a recommendation or license to use the product or any information that conflicts with any patent or intellectual property of MEGGLE or others and any such determination of use is the sole responsibility of the user. © MEGGLE