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Sieve data – spray-dried lactose			
	Lactose type	FlowLac® 90	FlowLac® 100
		specified/typical	specified/typical
Particle size distribution	< 32 µm	NMT 5%/2%	NMT 10%/5%
Method:	< 100 µm	25–40%/29%	20–45%/32%
Air jet sieving	< 200 µm	NLT 85%/91%	NLT 80%/87%
	< 250 µm	/99%	/97%

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
Figure 5: Sorption-desorption isotherm of spray-dried lactose, using FlowLac® 90 as an example.

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Functional related characteristics

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Packaging and shelf life
FlowLac®

	Size	Material	Shelf life
FlowLac® 90	25 kg	Corrugated carton box with an aluminium laminated inliner	36 months
FlowLac® 100			24 months
FlowLac® 100		Paper bag with PE-EVOH-PE inliner	18 months

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Literature

- [1] Meeus, L. (2011). Direct Compression versus Granulation. *Pharmaceutical Technology*, 23(3).
- [2] Kristensen, H. G., & Schaefer, T. (1987). Granulation: A Review on Pharmaceutical Wet-Granulation. *Drug Development and Industrial Pharmacy*, 13(4–5), 803–872.
- [3] Mîinea, L. A., Mehta, R., Kallam, M., Farina, J. A., & Deorkar, N. (2011). Evaluation and Characteristics of a New Direct Compression Performance Excipient, 35(3).
- [4] Gohel, M. C., & Jogani, P. D. (2005). A review of co-processed directly compressible excipients. *Journal of pharmacy & pharmaceutical sciences: a publication of the Canadian Society for Pharmaceutical Sciences, Société canadienne des sciences pharmaceutiques*, 8(1), 76–93.
- [5] Roos, Y. H. (2002). Importance of glass transition and water activity to spray drying and stability of dairy powders. *Le Lait*, 82(4), 475–484.
- [6] RUANGCHAYAJATUPORN, J., AMORNSAKCHAI, T., SINCHAIPANID, N., & MITREVEJ, A. (n.d.). Compaction behavior and optimization of spray-dried lactose with various amorphous content. *Journal of drug delivery science and technology*, 21(2), 175–181.

MEGGLE App:

MEGGLE Consultant

MEGGLE Group Wasserburg
BG Excipients & Technology
Meggelstrasse 6–12
83512 Wasserburg
Germany

Phone +49 8071 73 476
Fax +49 8071 73 320
service.pharma@meggle.de
www.meggle-pharma.com

MEGGLE warrants that its products conform to MEGGLE's written specification and makes no other expressed or implied warranties or representations. For any specific usage, the determination of suitability of use or application of MEGGLE products is the sole responsibility of the user. The determination of the use, application, and compliance of this product with regard to any national, regional, or local laws and/or regulations is the sole responsibility of the user, and MEGGLE makes no representation with regards to same. Nothing herein shall be construed as a recommendation or license to use the product or any information that conflicts with any patent or intellectual property of MEGGLE or others and any such determination of use is the sole responsibility of the user. © MEGGLE